


# “Solution missing” Error - tiamo

This help note describes how to correct ‘solution missing’ error in tiamo software.

*Technology: Titration*

## How do you know you have this problem?

Due to one of two possibilities, you have encountered the following error:



## What can I try to resolve it?

The way to find the source of the error is to check your Configuration > Titrants/Solutions:

Titrants/Solutions							
	Solution name ^	Concentration	Cylinder volume	Type	Dosing device	Titer	Date titer det.
1	acetate buffer	1 mol/L	50	IDU	859_2 / D2	1 mol/L	2013-08-15 09:2...
2	Al(NO3)	0.445 mol/L	20	IDU	859_2 / D1	0.445 mol/L	2013-08-15 09:1...
3	ammonium sulfate	1 mol/L	20	IDU		1 mol/L	2013-07-10 09:4...
4	Comp 5	1.0 mol/L	10	IDU		5.4 mol/L	2013-04-22 12:5...
5	hexadecylpyridini...	0.02 mol/L	20	IDU		0.02	2013-03-25 13:0...
6	Hydrochloric Acid	0.1 mol/L	2	IDU		5.02 mg/g	2012-05-30 14:5...
7	NaOH	1 mol/L	20	IDU		1 mol/L	2013-07-01 11:1...
8	Perchloric	0.1 mol/L	50	IDU		1	2013-05-31 10:5...
9	Perchloric Acid	0.1 mol/L	20	IDU	907_2 / D1	0.1	2013-07-15 12:2...
10	Silver Nitrate	0.141 mol/L	20	IDU	907_2 / D2	0.141	2013-07-08 12:0...
11	TSI	1 M	20	IDU		1 M	2013-01-10 13:2...
12	Water	1.000 mol/L	5	IDU		1.000	2013-02-14 16:2...

- IMPORTANT – On line 10, we have a solution named “Silver Nitrate”.
  - ‘Type’ in 4<sup>th</sup> column: IDU
    - This means solution is on Intelligent Dosing Unit, i.e. a dosing unit with data chip. This could also say IEU which stands for Intelligent Exchange Unit.

- 'Dosing device' in 5<sup>th</sup> column: 907\_2/D2
  - This means the dosing unit is currently connected to device 907\_2 on MSB port 2.

If your solution is being actively recognized and read, i.e. there is an entry in the 'Dosing device' column in Titrants/Solutions, continue with **Scenario 1**.


If your solution is NOT being actively recognized and read, i.e. there is NOT an entry in the 'Dosing device' column in Titrants/Solutions, continue with **Scenario 2**.

## Scenario 1

---

The cause of this error is a discrepancy between how your titrator is configured (physical connection of dosing device) and how the method is programmed.

1. Click Method tab
2. Click File > Open – select the method you're trying to run showing error
3. Double click on the command causing this error
  - a. In this example, the command name is 'DET U' per the error message on page 1



Command name: DET U

Potentiometric evaluation | Additional evaluations | Additional measured values

General/Hardware | Start conditions | Titration parameters | Stop conditions

**Device**

Device name: 907\_2

Device type: 907 Titrande

**Dosing device**

Dosing device: 1

Solution: Silver Nitrate

**Sensor**

Measuring input: 1

Sensor: Metal electrode

Temperature measurement: automatic

**Stirrer**

Stirrer: 4

Stirring rate: 6

Switch off automatically

OK Cancel


- IMPORTANT – Under ‘Dosing device’ section, Dosing device 1 is programmed for Solution ‘Silver Nitrate’. We know from Configuration > Titrants/Solutions, Silver Nitrate is on dosing device 2 (907\_2/D2).
- 4. Update the method – changing dosing device designation from 1 to 2
- 5. Click OK, then File > Save to save changes made to method
- 6. Alternatively, you can change the physical connection of the dosing device to the titrator – moving dosing device cable from MSB port 2 to MSB port 1  
(Note: Exit tiemo software before disconnecting dosing devices from titrator)

## Scenario 2

The cause of this error is the dosing unit or exchange unit is not properly communicating with the titrator; likely due to faulty data chip.

To overcome, we have two options:

1. 'Submit a request' for further assistance from Metrohm Technical Support at [support.metrohmusa.com](https://support.metrohmusa.com) – requesting assistance with procuring a new data chip
2. Program your method to run with Solution: 'not defined'
  - a. Click Method tab
  - b. Click File > Open – select the method you're trying to run showing error
  - c. Double click on the command causing this error
    - i. In this example, the command name is 'DET U' per the error message on page 1
  - d. Confirm physical connection of dosing device to rear of titrator (i.e. MSB port number connection)
  - e. Ensure 'Dosing device' number programmed in method matches connection to MSB port number on rear of titrator
  - f. Use drop down menu under 'Solution' to select 'not defined'



DET U - DET U

Command name: DET U

Potentiometric evaluation | Additional evaluations | Additional measured values

General/Hardware | Start conditions | Titration parameters | Stop conditions

**Device**

Device name: 809\_1

Device type: 809 Titrando

**Dosing device**

Dosing device: 1

Solution: not defined

**Sensor**

Measuring input: 1

Sensor: Metal electrode

Temperature measurement: automatic

**Stirrer**

Stirrer: 4

Stirring rate: 6

Switch off automatically

OK Cancel

- g. Click OK, then File > Save to save changes made to method

## Other Ideas

---

To remove defective data chip:

- The dosing unit chip resides on the top facing portion of the dosing unit



- Remove dosing unit housing from cylinder/distributor (press button above "Vent" and rotate counterclockwise)
- Using your finger, push up on underside of data chip from the inside of the housing

'Submit a request' for further assistance from Metrohm Technical Support at [support.metrohmusa.com](https://support.metrohmusa.com).